

# (8-001) Fehlermeldung bei Materialeingabe

## Fehlermeldung:

The screenshot shows a software interface with several input fields. The 'Material' field contains the text '21151 Filterkuchen mit 3% SiC Innogutschein Düro wdH'. A modal dialog box is displayed in the center, containing the message: 'Validation Failed for one or more entities. See 'EntityValidationErrors' property for more details.' The dialog has an 'OK' button at the bottom right. The background interface includes fields for 'Ordner d. Messung' (M2106170950), 'Gruppe (Ordner)' (P\_Max (G1901090646)), 'Messung' (M2106170950), 'Material' (21151 Filterkuchen mit 3% SiC Innogutschein Düro wdH), 'Probekörper', 'Anmerkung', 'Geräte-ID' (Erhitzungsmikrosko), and 'Datum' (6/17/2021 9:50:13).

“Validation failed for one or more entities. ...”

## Phänomen:

Nach der Eingabe der Materialbezeichnung wird die oben gezeigte Fehlermeldung angezeigt.

## Beobachtung(en), Annahmen:

- Die Materialbezeichnung ist ZU lang.
- Es sind hier maximal 50 Zeichen zugelassen.

## Abhilfe:

- Die Materialbezeichnung wird unter der Beachtung der “50-Zeichen-Regel” eingegeben.
- Weitere Erklärungen zum Material, der Probenentnahme, der Probenmaterialvorbereitung (z.B. Veraschungstemperatur, Mahlgrad, Partikelgröße etc.) werden im Feld ‘Anmerkung’ untergebracht.

## Weiteres

Fehler und Fehlermeldungen können (bei Versionen 3.1.6 und 3.2.0) ebenfalls auftreten, wenn eine zu lange Messungsbezeichnung eingegeben wird/wurde.

S. a. *„Messung nicht abgespeichert“*